

# CSF FLOW CONTROL VALVE BURR-HOLE

Sanitary Authorization Number 222E82 SSA



## ► PRESENTATION:

Individual. [Sterile and disposable product]

**YOUR WELLNESS IS OUR MOTIVATION**

## ► DESCRIPTION

This valve is indicated in the treatment of hydrocephalus in all its varieties, as there is increased intracranial pressure or when you want to improve brain function through a hypotension regime.

## ► THIS DEVICE CONSIST OF

A biconvex camera for adult and pediatric sizes.

A flat base for neonatal size camera.



A radiopaque ventricular catheter 15cm in length for adult and pediatric sizes.

A radiopaque ventricular catheter 20 cm in length for neonatal size.

A peritoneal catheter 90cm in length for adult and pediatric sizes.

A peritoneal catheter 75cm in length for neonatal size.

Three straight connectors 14mm in length.





# CSF FLOW CONTROL VALVE BURR-HOLE




Sanitary Authorization Number 222E82 SSA

YOUR WELLNESS IS OUR MOTIVATION

## BIOMED's CODE

CODE	SIZE	PRESSURE	OPERATING PRESSURE*	CODE
VDVP-A16	Adult	High	95-125	●●●
VDVP-M16	Adult	Medium	60-80	●●
VDVP-B16	Adult	Low	30-45	●
VDVP-A14	Pediatric	High	95-125	●●●
VDVP-M14	Pediatric	Medium	60-80	●●
VDVP-B14	Pediatric	Low	30-45	●
VDVP-A12	Pediatric	High	95-125	●●●
VDVP-M12	Pediatric	Medium	60-80	●●
VDVP-B12	Pediatric	Low	30-45	●
VDVP-NA	Neonatal	High	95-125	●●●
VDVP-NM	Neonatal	Medium	60-80	●●
VDVP-MB	Neonatal	Low	30-45	●

\*OPERATING PRESSURE IN: mm H<sub>2</sub>O



Manufactures in México by  
**Biomédica Mexicana, S.A. de C.V.**

Calzada Ermita Iztapalapa, No. 855  
Colonia Santa Isabel Industrial  
Delegación Iztapalapa, C.P. 09820  
Ciudad de México, Mexico

+52 (55) 5581 5919 / +52 (55) 5582 7036

ventas@biomedicamexicana.com.mx

@Biomedica\_mex

/BiomedicaMexicana

www.biomedicamexicana.com.mx

