

EMBOLECTOMY CATHETERS

Sanitary Authorization Number 55E83 SSA



► PRESENTATION:

Individual. [Sterile and disposable product]

YOUR WELLNESS IS OUR MOTIVATION

► DESCRIPTION

It is indicated for removing clots formed in the arteries. It has referenced lines each 10cm from its distal end.

► THIS DEVICE CONSIST OF

Unilumen tube of polyvinylchloride.

A latec balloon for removal of clot.

A translucent luer adapter at its proximal end.

A stainless steel stylet running along the lumen.

BIOMED'S CODE

CODE	LENGHT (cm)	FRENCH	COLOR	MAXIMUM CAPACITY OF LIQUID(ml)
SOEM-403	40	3	●	0.20
SOEM-404	40	4	●	0.20
SOEM-405	40	5	○	0.20
SOEM-406	40	6	●	0.20
SOEM-407	40	7	●	0.20
SOEM-803	80	3	●	0.20
SOEM-804	80	4	●	0.20
SOEM-805	80	5	○	0.20
SOEM-806	80	6	●	0.20
SOEM-807	80	7	●	0.20



Manufactures in México by
Biomédica Mexicana, S.A. de C.V.

Calzada Ermita Iztapalapa, No. 855
Colonia Santa Isabel Industrial
Delegación Iztapalapa. C.P. 09820
Ciudad de México, Mexico

+52 (55) 5581 5919 / +52 (55) 5582 7036

ventas@biomedicamexicana.com.mx

@Biomedica_mex

/BiomedicaMexicana

www.biomedicamexicana.com.mx