

# EXTERNAL DRAINAGE SYSTEM

Sanitary Authorization Number 2724C2012 SSA



## PRESENTATION:

Individual. [Sterile and disposable product]

YOUR WELLNESS IS OUR MOTIVATION

## BIOMED'S CODE





CODE	MONITORING
EDVE-CM	YES
EDVE-SM	NO

## DESCRIPTION

This equipment is intended to reduce and monitor intracranial pressure and cerebrospinal fluid, also works as a temporary drainage in patients with CSF shunt valves for presenting infection.

## THIS DEVICE CONSIST OF

	EDVE-CM	EDVE-SM
A ventricular catheter radiopaque 30cm long, with suture ring and stainless steel stylet.	<input type="radio"/>	<input type="radio"/>
Luer integrated system and ventricular drainage catheter	<input type="radio"/>	<input type="radio"/>
A collection bag to 700ml graduated with antibacterial and hydrophobic filter and security keys.	<input type="radio"/>	<input type="radio"/>
A drip chamber graduated to 70ml with key 4-way output, as well as antibacterial and hydrophobic filter.	<input type="radio"/>	<input type="radio"/>
An anti-reflux valve.	<input type="radio"/>	<input type="radio"/>
4 way stopcock in the patient's piping.	<input type="radio"/>	<input type="radio"/>
Stainless Steel curved trocar.	<input type="radio"/>	<input type="radio"/>
170cm of PVC tubing.	<input type="radio"/>	<input type="radio"/>
Scale graduated in mm Hg and in cm H <sub>2</sub> O.	<input type="radio"/>	<input type="radio"/>
Needleless injection site.	<input type="radio"/>	<input type="radio"/>
Flat stopcocks in the patient's piping.	<input type="radio"/>	<input type="radio"/>



Manufactures in México by **Biomédica Mexicana, S.A. de C.V.**

Calzada Ermita Iztapalapa, No. 855  
 Colonia Santa Isabel Industrial  
 Delegación Iztapalapa. C.P. 09820  
 Ciudad de México, Mexico

+52 (55) 5581 5919 / +52 (55) 5582 7036

ventas@biomedicamexicana.com.mx

@Biomédica\_mex

/BiomédicaMexicana

www.biomedicamexicana.com.mx