

LACRIMAL INTUBATION EQUIPMENT

Sanitary Authorization Number 0967C2016 SSA



► PRESENTATION:

Individual. [Sterile and disposable product]

YOUR WELLNESS IS OUR MOTIVATION

BIOMED's CODE

CODE
EILA




► DESCRIPTION

Indicated for the treatment of obstruction of the tear ducts. **BIOMED's lacrimal intubation equipment** is a device indicated to prevent the consequences of this obstruction, such a chronic infections, eliminating the impediment of normal drainage and restoring it.

► THIS DEVICE CONSIST OF

Two stainless steel probes 0.60mm in diameter and 17cm in length.

A silicone tube 30cm long.



Manufactures in México by
Biomédica Mexicana, S.A. de C.V.

Calzada Ermita Iztapalapa, No. 855
Colonia Santa Isabel Industrial
Delegación Iztapalapa. C.P. 09820
Ciudad de México, Mexico

+52 (55) 5581 5919 / +52 (55) 5582 7036

ventas@biomedicamexicana.com.mx

@Biomédica_mex

/BiomédicaMexicana

www.biomedicamexicana.com.mx