

OMMAYA'S CSF RESERVOIR

Sanitary Authorization Number 164E84 SSA



► PRESENTATION:

Individual. [Sterile and disposable product]

YOUR WELLNESS IS OUR MOTIVATION

BIOMED's CODE

CODE	SIZE
RLCR-P	PEDIATRIC 12mm
RLCR-A	ADULT 24mm

► DESCRIPTION

Suitable when you want to access the cerebral ventricles through simple maneuvers puncture in the dome, which can be taken CSF samples, temporarily decompress the ventricular system and administer drugs or contrast media.

► THIS DEVICE CONSIST OF




A dome-shaped reservoir with a flat base.

A ventricular catheter radiopaque 15cm in length.

A polypropylene elbow.

A straight connector 14mm in length.

A stainless steel stylet.



Manufactures in México by
Biomédica Mexicana, S.A. de C.V.

Calzada Ermita Iztapalapa, No. 855
Colonia Santa Isabel Industrial
Delegación Iztapalapa. C.P. 09820
Ciudad de México, Mexico

+52 (55) 5581 5919 / +52 (55) 5582 7036

ventas@biomedicamexicana.com.mx

@Biomedica_mex

/BiomedicaMexicana

www.biomedicamexicana.com.mx