

# UNIVERSAL SUBDURAL DRAINAGE BAG

Sanitary Authorization Number 1128C2009 SSA



## ► PRESENTATION:

Individual. [Sterile and disposable product]

YOUR WELLNESS IS OUR MOTIVATION

## BIOMED's CODE

CODE  
BDSU

## ► DESCRIPTION

The objective of this device is to drain cumulative flow product of a hematoma, is designed to be placed through trephine holes or by craniotomy. All the device are 100% silicone.

## ► THIS DEVICE CONSIST OF


Bag or reservoir with a capacity of 250ml.

Subdural catheter 14Fr and 30cm long; rounded tip in the first 5cm in intracranial extreme.

Suture ring for attachment to adjacent tissues.

Connector polypropylene.

Drain connector in the bag outlet valve for emptying or key.



Manufactures in México by  
**Biomédica Mexicana, S.A. de C.V.**

Calzada Ermita Iztapalapa, No. 855  
Colonia Santa Isabel Industrial  
Delegación Iztapalapa. C.P. 09820  
Ciudad de México, Mexico

+52 (55) 5581 5919 / +52 (55) 5582 7036

ventas@biomedicamexicana.com.mx

@Biomedica\_mex

/BiomedicaMexicana

www.biomedicamexicana.com.mx